

Product datasheet for **CD563179**

Tissue Genomic DNA, Tonsil

Product data:

Product Type:	Genomic DNA
Disease State:	Other Disease
Diagnosis Category:	Hyperplasia
Tissue:	Tonsil
Frozen Sample ID:	FR00000DC6
Case ID:	CU0000000387
Prep:	A
Vial ID:	DN000000BB
Tissue of Origin:	Tonsil
Site of Finding:	Tonsil
Appearance:	Lesional
Sample Diagnosis (from Pathology Verification):	Follicular hyperplasia of tonsil
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	10% squamous epithelium, 45% reactive secondary follicles, 10% interfollicular lymphoid cells, 35% fibrous tissue
Bioanalyzer Ratio (28S/18S):	2.05



[View online »](#)

CASE Diagnosis (from medical center path report):

Follicular hyperplasia of tonsil


Age: 31

Gender: Female


Storage: Store the total DNA -80°C.

Stability: These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Product images:



4x Image



20x Image