

Product datasheet for **CD563180**

Tissue Genomic DNA, Prostate

Product data:

| | |
|---|--|
| Product Type: | Genomic DNA |
| Disease State: | Other Disease |
| Diagnosis Category: | Benign Prostatic Hyperplasia (BPH) |
| Tissue: | Prostate |
| Frozen Sample ID: | FR0001AC12 |
| Case ID: | CI0000006008 |
| Prep: | A |
| Vial ID: | DN00001418 |
| Tissue of Origin: | Prostate |
| Site of Finding: | Prostate |
| Appearance: | Lesional |
| Sample Diagnosis (from Pathology Verification): | Hyperplasia of prostate |
| Normal: | 0% |
| Lesional: | 100% |
| Tumor: | 0% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 0% |
| Necrosis: | 0% |
| Pathology Verification notes from H&E review: | Non Tumor Structures: 25% Epithelium, 75% Fibromuscular stroma |
| Bioanalyzer Ratio (28S/18S): | 2 |



[View online »](#)

| | |
|--|--|
| CASE Diagnosis (from medical center path report): | Hyperplasia of prostate |
| Age: | 59 |
| Gender: | Male |
| Race: | Not reported |
| Storage: | Store the total DNA -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |