

Product datasheet for **CD563313**

Tissue Genomic DNA, Liver

Product data:

Product Type:	Genomic DNA
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Liver
Frozen Sample ID:	FR0002E1E9
Case ID:	CI0000010635
Prep:	B
Vial ID:	DN00000BD5
Tissue of Origin:	Liver
Site of Finding:	Liver
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Inflammation: Mild Lymphoplasmacytic infiltrate. Non Tumor Structures: 85% Lobules, 15% Triads
Bioanalyzer Ratio (28S/18S):	2


[View online »](#)

CASE Diagnosis (from medical center path report): Adenocarcinoma of colon, metastatic

Age: 51

Gender: Female

Race: White or Caucasian

Tumor Grade: Not applicable

Minimum Stage Grouping: Not applicable

T (Extent of Primary Tumor): Not applicable


N (Lymph node metastasis): Not applicable

M (Distant metastasis): Not applicable


Storage: Store the total DNA -80°C.

Stability: These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Product images:



DNA - 4x.



DNA - 20x.



DNA PCR Image