

Product datasheet for **CD563945**

Tissue Genomic DNA, Kidney

Product data:

| | |
|---|---|
| Product Type: | Genomic DNA |
| Disease State: | Cancer |
| Diagnosis Category: | Renal Cell Carcinoma |
| Tissue: | Kidney |
| Frozen Sample ID: | FR00001FE7 |
| Case ID: | CI0000000038 |
| Prep: | B |
| Vial ID: | DN000000D7 |
| Tissue of Origin: | Kidney |
| Site of Finding: | Kidney |
| Appearance: | Tumor |
| Sample Diagnosis (from Pathology Verification): | Carcinoma of kidney, renal cell, clear cell |
| Normal: | 0% |
| Lesional: | 0% |
| Tumor: | 90% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 10% |
| Necrosis: | 0% |
| Bioanalyzer Ratio (28S/18S): | 1.95 |
| CASE Diagnosis (from medical center path report): | Carcinoma of kidney, renal cell, clear cell |
| Age: | 45 |



[View online »](#)

| | |
|-------------------------------------|--|
| Gender: | Female |
| Race: | Black or African American |
| Tumor Grade: | Fuhrman G1: Nuclei round, uniform |
| Minimum Stage Grouping: | I |
| T (Extent of Primary Tumor): | T1 |
| N (Lymph node metastasis): | NX |
| M (Distant metastasis): | MX |
| Storage: | Store the total DNA -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |