

Product datasheet for **CD564108**

Tissue Genomic DNA, Lung

Product data:

| | |
|---|----------------------------------|
| Product Type: | Genomic DNA |
| Disease State: | Cancer |
| Diagnosis Category: | Lung Cancer |
| Tissue: | Lung |
| Frozen Sample ID: | FR00022D3B |
| Case ID: | CI0000009032 |
| Prep: | A |
| Vial ID: | DN000002A3 |
| Tissue of Origin: | Lung |
| Site of Finding: | Lung |
| Appearance: | Tumor |
| Sample Diagnosis (from Pathology Verification): | Carcinoma of lung, squamous cell |
| Normal: | 0% |
| Lesional: | 0% |
| Tumor: | 70% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 5% |
| Necrosis: | 25% |
| Bioanalyzer Ratio (28S/18S): | 1.97 |
| CASE Diagnosis (from medical center path report): | Carcinoma of lung, squamous cell |
| Age: | 58 |



[View online »](#)

| | |
|-------------------------------------|--|
| Gender: | Male |
| Race: | Not reported |
| Tumor Grade: | AJCC G3: Poorly differentiated |
| Minimum Stage Grouping: | IA |
| T (Extent of Primary Tumor): | T1 |
| N (Lymph node metastasis): | N0 |
| M (Distant metastasis): | MX |
| Diagnostic Test Results: | Mucicarmine,Negative |
| Storage: | Store the total DNA -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |