

Product datasheet for **CD564580**

Tissue Genomic DNA, Colon

Product data:


| | |
|---|--|
| Product Type: | Genomic DNA |
| Disease State: | Normal |
| Diagnosis Category: | Normal |
| Tissue: | Colon |
| Frozen Sample ID: | FR0003DF71 |
| Case ID: | CU0000014093 |
| Prep: | B |
| Vial ID: | DN00000E57 |
| Tissue of Origin: | Colon |
| Site of Finding: | Colon |
| Appearance: | Normal |
| Sample Diagnosis (from Pathology Verification): | Within normal limits |
| Normal: | 100% |
| Lesional: | 0% |
| Tumor: | 0% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 0% |
| Necrosis: | 0% |
| Pathology Verification notes from H&E review: | Non Tumor Structures: 75% Mucosa, 20% Submucosa, 5% Muscle |
| Bioanalyzer Ratio (28S/18S): | 2.03 |




[View online »](#)

| | |
|---|--|
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of colon |
| Age: | 81 |
| Gender: | Male |
| Tumor Grade: | Not applicable |
| Minimum Stage Grouping: | Not applicable |
| T (Extent of Primary Tumor): | Not applicable |
| N (Lymph node metastasis): | Not applicable |
| M (Distant metastasis): | Not applicable |
| Storage: | Store the total DNA -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |

Product images:



DNA - 4x.



DNA - 20x.



DNA PCR Image