

## Product datasheet for **CD564803**

### Tissue Genomic DNA, Liver

#### Product data:

Product Type:	Genomic DNA
Disease State:	Other Disease
Diagnosis Category:	Cirrhosis
Tissue:	Liver
Frozen Sample ID:	FR00035B4E
Case ID:	CU0000011777
Prep:	C
Vial ID:	DN000006BE
Tissue of Origin:	Liver
Site of Finding:	Liver
Appearance:	Lesional
Sample Diagnosis (from Pathology Verification):	Cirrhosis of liver
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Inflammation: Mild Mixed inflammatory infiltrate. Non Tumor Structures: 95% Lobules, 5% Triads. Other Features/Comments: Chronic active hepatitis with cirrhosis, compatible with Hepatitis B infection (per path report)
Bioanalyzer Ratio (28S/18S):	1.72



[View online »](#)

<b>CASE Diagnosis (from medical center path report):</b>	Carcinoma of liver, hepatocellular
<b>Age:</b>	58
<b>Gender:</b>	Female
<b>Tumor Grade:</b>	Not Reported
<b>Minimum Stage Grouping:</b>	I
<b>T (Extent of Primary Tumor):</b>	T1
<b>N (Lymph node metastasis):</b>	NX
<b>M (Distant metastasis):</b>	MX
<b>Storage:</b>	Store the total DNA -80°C.
<b>Stability:</b>	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.