

## Product datasheet for **CD564843**

### Tissue Genomic DNA, Lung

#### Product data:

|   |                                  |
|---|----------------------------------|
| Product Type:                                     | Genomic DNA                      |
| Disease State:                                    | Cancer                           |
| Diagnosis Category:                               | Lung Cancer                      |
| Tissue:   | Lung                             |
| Frozen Sample ID:                                 | FR000015E4                       |
| Case ID:  | CU0000000597                     |
| Prep:   | B                                |
| Vial ID:  | DN00000304                       |
| Tissue of Origin:                                 | Lung                             |
| Site of Finding:                                  | Lung                             |
| Appearance:                                       | Tumor                            |
| Sample Diagnosis (from Pathology Verification):   | Carcinoma of lung, squamous cell |
| Normal:   | 0%                               |
| Lesional:   | 0%                               |
| Tumor:  | 60%                              |
| Tumor Hypo/Acellular Stroma:                      | 0%                               |
| Tumor Hypercellular Stroma:                       | 15%                              |
| Necrosis:   | 25%                              |
| Bioanalyzer Ratio (28S/18S):                      | 1.92                             |
| CASE Diagnosis (from medical center path report): | Carcinoma of lung, squamous cell |
| Age:  | 71                               |



[View online »](#)

|                                     |  |
|-------------------------------------|--|
| <b>Gender:</b>                      | Male   |
| <b>Tumor Grade:</b>                 | AJCC G2: Moderately differentiated   |
| <b>Minimum Stage Grouping:</b>      | IIB  |
| <b>T (Extent of Primary Tumor):</b> | T3   |
| <b>N (Lymph node metastasis):</b>   | N0   |
| <b>M (Distant metastasis):</b>      | MX   |
| <b>Storage:</b>                     | Store the total DNA -80°C.   |
| <b>Stability:</b>                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |