

Product datasheet for **CF180016**

Carrier-free (BSA/glycerol-free) Mouse monoclonal E2-Crimson antibody, Clone OT13E12


Product data:

Product Type:	Tag Antibodies
Product Name:	Carrier-free (BSA/glycerol-free) Mouse monoclonal E2-Crimson antibody, Clone OT13E12
Clone Name:	OT13E12
Applications:	WB
Recommended Dilution:	WB: 1:500
Reactivity:	E2-Crimson
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	A E2-Crimson tagged recombinant protein expressed in HEK293T cell.
Isotype:	IgG1
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	26.2 kDa



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY E2-Crimson (JA1132-C07, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-E2-Crimson (1:500).