

Product datasheet for CF180148

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Carrier-free (BSA/glycerol-free) Mouse monoclonal tdTomato antibody, Clone OTI4B3

Product data:

Product Type: Tag Antibodies

Product Name: Carrier-free (BSA/glycerol-free) Mouse monoclonal tdTomato antibody, Clone OTI4B3

Clone Name: OTI4B3

Applications: WB

Recommended Dilution:WB: 1:500Reactivity:tdTomatoHost:Mouse

Clonality: Monoclonal

Immunogen: A tdTomato tagged recombinant protein expressed in E.coli.

Isotype: IgG1

Formulation: Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)

Reconstitution Method: For reconstitution, we recommend adding 100uL distilled water to a final antibody

concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)

Conjugation: Unconjugated Storage: Store at -20°C.


Stability: Stable for 1 year at -20°C from date of shipment.

Predicted Protein Size: 54.2 kDa





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY tdTomato (JB1016_D07, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-tdTomato (1:500).