

Product datasheet for **CF500878**

PASK Mouse Monoclonal Antibody [Clone ID: OTI11D8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI11D8
Applications:	IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human PASK (NP_055963) produced in HEK293T cell.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	142.9 kDa
Gene Name:	PAS domain containing serine/threonine kinase
Database Link:	NP_055963 Entrez Gene 23178 Human Q96RG2
Background:	PAS domains regulate the function of many intracellular signaling pathways in response to both extrinsic and intrinsic stimuli. PASK is an evolutionarily conserved protein present in yeast, flies, and mammals.




[View online »](#)


Synonyms: PASKIN; STK37

Protein Families: Druggable Genome, Protein Kinase, Stem cell - Pluripotency

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PASK ([RC209096], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PASK. Positive lysates [LY414744] (100ug) and [LC414744] (20ug) can be purchased separately from OriGene.



Anti-PASK mouse monoclonal antibody ([TA500878]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PASK ([RC209096]).