

Product datasheet for CF501845

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MICAL1 Mouse Monoclonal Antibody [Clone ID: OTI6H10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6H10

Applications: FC, IHC, WB

Recommended Dilution: WB 1:500~2000, FLOW 1:100, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MICAL1 (NP_073602) produced in HEK293T

cell.

Formulation: Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)

Reconstitution Method: For reconstitution, we recommend adding 100uL distilled water to a final antibody

concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 117.7 kDa

Gene Name: microtubule associated monooxygenase, calponin and LIM domain containing 1

Database Link: NP 073602

Entrez Gene 64780 Human


Q8TDZ2

Synonyms: MICAL; MICAL-1; NICAL







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MICAL1 ([RC208308], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MICAL1. Positive lysates [LY402940] (100ug) and [LC402940] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-MICAL1 monoclonal antibody.




Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-MICAL1 mouse monoclonal antibody ([TA501845]) at 1:150 dilution.






Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-MICAL1 mouse monoclonal antibody ([TA501845]) at 1:150 dilution.



HEK293T cells transfected with either [RC208308] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-MICAL1 antibody ([TA501845]), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-MICAL1 antibody ([TA501845]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).