

## Product datasheet for **CF502684**

### **RNF156 (MGRN1) Mouse Monoclonal Antibody [Clone ID: OTI3C3]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI3C3
<b>Applications:</b>	FC, IF, WB
<b>Recommended Dilution:</b>	WB 1:500~2000, IF 1:100, FLOW 1:100
<b>Reactivity:</b>	Human, Monkey, Mouse, Rat
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human MGRN1 (NP_056061) produced in HEK293T cell.
<b>Formulation:</b>	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
<b>Reconstitution Method:</b>	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Predicted Protein Size:</b>	63 kDa
<b>Gene Name:</b>	mahogunin ring finger 1
<b>Database Link:</b>	<a href="#">NP_056061</a> <a href="#">Entrez Gene 709272 Monkey</a> <a href="#">Entrez Gene 23295 Human</a> <a href="#">O60291</a>
<b>Background:</b>	Mahogunin (MGRN1) is a C3HC4 RING-containing protein with E3 ubiquitin ligase activity in vitro. [supplied by OMIM]




[View online »](#)


Synonyms: RNF156

Protein Pathways: Ubiquitin mediated proteolysis


**Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY MGRN1 (Cat# [RC208284], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MGRN1(Cat# [TA502684]). Positive lysates [LY414714] (100ug) and [LC414714] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-MGRN1 monoclonal antibody.




Anti-MGRN1 mouse monoclonal antibody ([TA502684]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MGRN1 ([RC208284]).



HEK293T cells transfected with either [RC208284] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-MGRN1 antibody ([TA502684]), and then analyzed by flow cytometry.



Flow cytometric Analysis of HeLa cells, using anti-MGRN1 antibody ([TA502684]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-MGRN1 antibody ([TA502684]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).