

## Product datasheet for **CF502914**

### NEU2 Mouse Monoclonal Antibody [Clone ID: OTI5E1]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI5E1
Applications:	IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NEU2 (NP_005374) produced in HEK293T cell.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42.1 kDa
Gene Name:	Homo sapiens neuraminidase 2 (NEU2), mRNA.
Database Link:	<a href="#">NP_005374 Entrez Gene 4759 Human</a>
Background:	This gene belongs to a family of glycohydrolytic enzymes which remove sialic acid residues from glycoproteins and glycolipids. Expression studies in COS7 cells confirmed that this gene encodes a functional sialidase. Its cytosolic localization was demonstrated by cell fractionation experiments. [provided by RefSeq]
Synonyms:	SIAL2




[View online »](#)

Protein Pathways: Other glycan degradation, Sphingolipid metabolism

### Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NEU2 ([RC219858], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NEU2. Positive lysates [LY401652] (100ug) and [LC401652] (20ug) can be purchased separately from OriGene.



Anti-NEU2 mouse monoclonal antibody ([TA502914]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY NEU2 ([RC219858]).