

## Product datasheet for **CF503496**

### DTYMK Mouse Monoclonal Antibody [Clone ID: OT11C4]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OT11C4
Applications:	FC, IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human DTYMK(NP_036277) produced in HEK283T cell.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	23.6 kDa
Gene Name:	deoxythymidylate kinase
Database Link:	<a href="#">NP_036277</a> <a href="#">Entrez Gene 1841 Human P23919</a>
Synonyms:	CDC8; PP3731; TMPK; TYMK
Protein Families:	Druggable Genome




[View online »](#)

Protein Pathways: Metabolic pathways, Pyrimidine metabolism


### Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DTYMK ([RC201228], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DTYMK. Positive lysates [LY415933] (100ug) and [LC415933] (20ug) can be purchased separately from OriGene.



Anti-DTYMK mouse monoclonal antibody ([TA503496]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY DTYMK ([RC201228]).



HEK293T cells transfected with either [RC201228] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DTYMK antibody ([TA503496]), and then analyzed by flow cytometry.