

Product datasheet for **CF504582**

CYB5R1 Mouse Monoclonal Antibody [Clone ID: OTI5G9]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5G9
Applications:	IF, WB
Recommended Dilution:	WB 1:1000, IF 1:100
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CYB5R1(NP_057327) produced in HEK293T cell.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	33.9 kDa
Gene Name:	cytochrome b5 reductase 1
Database Link:	NP_057327 Entrez Gene 72017 Mouse Entrez Gene 304805 Rat Entrez Gene 51706 Human Q9UHQ9
Synonyms:	B5R.1; B5R1; B5R2; humb5R2; NQO3A2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: Amino sugar and nucleotide sugar metabolism

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CYB5R1 ([RC205833], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CYB5R1. Positive lysates [LY402525] (100ug) and [LC402525] (20ug) can be purchased separately from OriGene.



Anti-CYB5R1 mouse monoclonal antibody ([TA504582]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY CYB5R1 ([RC205833]).