

Product datasheet for **CF507127**


Angiomotin (AMOT) Mouse Monoclonal Antibody [Clone ID: OTI1A8]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1A8
Applications:	IF, WB
Recommended Dilution:	WB 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human AMOT(NP_573572) produced in HEK293T cell.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	72.4 kDa
Gene Name:	angiomotin
Database Link:	NP_573572 Entrez Gene 27494 Mouse Entrez Gene 300289 Rat Entrez Gene 154796 Human Q4VCS5
Synonyms:	KIAA1071



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AMOT ([RC220347], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AMOT. Positive lysates [LY408854] (100ug) and [LC408854] (20ug) can be purchased separately from OriGene.



Anti-AMOT mouse monoclonal antibody ([TA507127]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY AMOT ([RC220347]).