

## Product datasheet for **CF800025**


### DPPA4 Mouse Monoclonal Antibody [Clone ID: OT1C1]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OT1C1
Applications:	IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 9-254 of human DPPA4 (NP_060659) produced in E.coli.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	33.4 kDa
Gene Name:	developmental pluripotency associated 4
Database Link:	<a href="#">NP_060659</a> <a href="#">Entrez Gene 55211 Human</a> <a href="#">Q7L190</a>
Synonyms:	2410091M23Rik



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DPPA4 ([RC205999], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DPPA4. Positive lysates [LY413234] (100ug) and [LC413234] (20ug) can be purchased separately from OriGene.



Anti-DPPA4 mouse monoclonal antibody ([TA800025]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY DPPA4 ([RC205999]).