

## Product datasheet for **CF800605**


### ZNF307 (ZKSCAN4) Mouse Monoclonal Antibody [Clone ID: OTI6G6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI6G6
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-316 of human ZKSCAN4 (NP_061983) produced in E.coli.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	61.4 kDa
Gene Name:	zinc finger with KRAB and SCAN domains 4
Database Link:	<a href="#">NP_061983</a> <a href="#">Entrez Gene 387032 Human</a> <a href="#">Q969J2</a>
Synonyms:	P1P373C6; ZNF307; ZNF427; ZSCAN36
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ZKSCAN4 (Cat# [RC204276], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZKSCAN4 (Cat# [TA800605]). Positive lysates [LY402737] (100ug) and [LC402737] (20ug) can be purchased separately from OriGene.