

Product datasheet for CF802325

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NDUFS2 Mouse Monoclonal Antibody [Clone ID: OTI2G2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2G2
Applications: IHC, WB
Recommended Dilution: WB 1:2000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 205-463 of human

NDUFS2 (NP_004541) produced in E.coli.

Formulation: Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)

Reconstitution Method: For reconstitution, we recommend adding 100uL distilled water to a final antibody

concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 49.1 kDa

Gene Name: NADH:ubiquinone oxidoreductase core subunit S2

Database Link: NP 004541

Entrez Gene 226646 MouseEntrez Gene 289218 RatEntrez Gene 4720 Human

075306





Background:


The protein encoded by this gene is a core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (complex I). Mammalian mitochondrial complex I is composed of at least 43 different subunits, 7 of which are encoded by the mitochondrial genome, and the rest are the products of nuclear genes. The iron-sulfur protein fraction of complex I is made up of 7 subunits, including this gene product. Complex I catalyzes the NADH oxidation with concomitant ubiquinone reduction and proton ejection out of the mitochondria. Mutations in this gene are associated with mitochondrial complex I deficiency. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]


Synonyms: CI-49


Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation,

Parkinson's disease

Product images:









Western blot analysis of A549 cell lysate (35ug) by using anti-NDUFS2 monoclonal antibody. Dilution: 1:500






Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-NDUFS2 mouse monoclonal antibody. ([TA802325]) Dilution: 1:150




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-NDUFS2 mouse monoclonal antibody. ([TA802325]) Dilution: 1:150




Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-NDUFS2 mouse monoclonal antibody. ([TA802325]) Dilution: 1:150






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-NDUFS2 mouse monoclonal antibody. ([TA802325]) Dilution: 1:150



Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-NDUFS2 mouse monoclonal antibody. ([TA802325]) Dilution: 1:150



Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-NDUFS2 mouse monoclonal antibody. ([TA802325]) Dilution: 1:150