


Product datasheet for **CF809207**

UEVLD Mouse Monoclonal Antibody [Clone ID: OTI11E2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI11E2
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 131-379 of human UEVLD(NP_060784) produced in E.coli.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42.1 kDa
Gene Name:	UEV and lactate/malate dehydrogenase domains
Database Link:	NP_060784 Entrez Gene 55293 Human Q8IX04
Synonyms:	ATTP; UEV3

[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY UEVLD (Cat# [RC215907], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UEVLD (Cat# [TA809207])(1:2000). Positive lysates [LY413157] (100ug) and [LC413157] (20ug) can be purchased separately from OriGene.