

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for CF811327

RNA5-8SN2 Mouse Monoclonal Antibody [Clone ID: OTI1C8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1C8
Applications:	IP
Recommended Dilution:	IP 1:30
Reactivity:	Streptococcus Pyogenes
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-1166 of human CAS9 produced in E.coli.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Formulation: Reconstitution Method:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose) For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific) Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography
Reconstitution Method: Purification:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific) Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Reconstitution Method: Purification: Conjugation:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific) Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) Unconjugated
Reconstitution Method: Purification: Conjugation: Storage:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific) Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) Unconjugated Store at -20°C as received.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Immunoprecipitation (IP) of Cas9 and Cas9nickase by using mouse monoclonal anti-CAS9 antibodies [TA811194] and [TA811327]. Mouse IgG control serves as the negative control. 293T cells were transfected with flag-tagged Cas9 overexpression plasmid, pCas-Guide (G) and pCas-Guide-nickase (N). Overexpression cell lysates were first precleared with ProteinG beads then incubated with beads crosslinked with anti-CAS9 antibody for overnight. The beads were then rinced with buffer and went through Western Blot analysis using anti-flag antibody ([TA180144]) (1:30).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US