

Product datasheet for **CF811365**

TXNRD1 Mouse Monoclonal Antibody [Clone ID: OTI5G3]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5G3
Applications:	WB
Recommended Dilution:	WB 1:500~2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TXNRD1 (NP_003321) produced in E.coli.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	Homo sapiens thioredoxin reductase 1 (TXNRD1), transcript variant 4, mRNA.
Database Link:	NP_003321 Entrez Gene 50493 MouseEntrez Gene 58819 RatEntrez Gene 7296 Human
Synonyms:	GRIM-12; TR; TR1; TRXR1; TXNR
Protein Families:	Druggable Genome
Protein Pathways:	Pyrimidine metabolism




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TXNRD1 ([RC202040], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TXNRD1 (1:2000).



Western blot analysis of extracts (35ug) from 3 different cell lines and 2 different tissue lysates by using anti-TXNRD1 monoclonal antibody (1:500).