

Product datasheet for **CF812076**


BCAS2 Mouse Monoclonal Antibody [Clone ID: OTI10H7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI10H7
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human BCAS2 (NP_005863) produced in E.coli.
Formulation:	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Reconstitution Method:	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	26 kDa
Gene Name:	BCAS2 pre-mRNA processing factor
Database Link:	NP_005863 Entrez Gene 68183 MouseEntrez Gene 295334 RatEntrez Gene 10286 Human O75934
Synonyms:	DAM1; Snt309; SPF27
Protein Pathways:	Spliceosome



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY BCAS2 (Cat# [RC205615], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BCAS2 (Cat# [TA812076]). Positive lysates [LY417014] (100ug) and [LC417014] (20ug) can be purchased separately from OriGene.