

# **Product datasheet for CF814079**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### PGP9.5 (UCHL1) Mouse Monoclonal Antibody [Clone ID: OTI1G7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1G7
Applications: ELISA

Recommended Dilution: ELISA 1:5000-10000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human UCHL1 (NP\_004172) produced in E.coli.

Formulation: Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)

**Reconstitution Method:** For reconstitution, we recommend adding 100uL distilled water to a final antibody

concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Conjugation:** Unconjugated

Predicted Protein Size: 24.8 kDa

Gene Name: ubiquitin C-terminal hydrolase L1

Database Link: NP 004172

Entrez Gene 7345 Human

P09936

**Background:** The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol

protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease. [provided by RefSeq, Sep

2009]

Synonyms: HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1





## PGP9.5 (UCHL1) Mouse Monoclonal Antibody [Clone ID: OTI1G7] - CF814079

**Protein Families:** Druggable Genome, Protease

**Protein Pathways:** Parkinson's disease