

Product datasheet for CF814434

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SARS-CoV-2 N Protein Mouse Monoclonal Antibody [Clone ID: OTI3E1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3E1

Applications: ELISA, WB

Recommended Dilution: WB 1:500-1000, ELISA 1:5000-10000

Reactivity: SARS-CoV-2

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Recombinant full length SARS-CoV-2 N protein produced in E.coli.

Formulation: Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)

Reconstitution Method: For reconstitution, we recommend adding 100uL distilled water to a final antibody

concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. (OriGene recommends Zeba Spin Desalting Columns, 7KMWCO from Thermo Scientific)

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.


Stability: Stable for 12 months from date of receipt.

Database Link: Entrez Gene 43740575 SARS-CoV-2






Product images:



Western blot analysis of recombinant full length SARS-CoV-2 N protein (Cat# [TP790189], 0.02ug) by using anti-SARS-CoV-2 N protein antibody (Cat# [TA814434]). (1:1000)



Western blot analysis of recombinant full length SARS-CoV-2 N protein (Cat# [TP750180], 0.02ug) by using anti-SARS-CoV-2 N protein antibody (Cat# [TA814434]). (1:1000)