


Product datasheet for CH14101-100

S100a10 Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Rat
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	S100 calcium binding protein A10
Database Link:	Entrez Gene 81778 Rat
Synonyms:	S100-A10, ANX2LG, CAL1L, CLP11, Calpactin-1 light chain


Product images:



P11 immunofluorescence in rat choroid plexus and FITC-conjugated donkey anti-chicken secondary antibodies.



[View online »](#)



P11 immunofluorescence in rat choroid plexus and FITC-conjugated donkey anti-chicken secondary antibodies.