


Product datasheet for CH23001-100

Nestin (NES) Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	nestin
Database Link:	Entrez Gene 10763 Human
Synonyms:	NES, Nbla00170


Product images:



A tissue section through an e13 mouse brain showing Nestin (green staining) in neural stem cells of the ventricular zone of the CNS. DAPI (blue staining) allows visualization of nuclei.



[View online »](#)



A tissue section through an e13 mouse brain showing Nestin (green staining) in neural stem cells of the ventricular zone of the CNS. DAPI (blue staining) allows visualization of nuclei.