


## Product datasheet for CH23004-200

### GAP43 Chicken Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	growth associated protein 43
Database Link:	<a href="#">Entrez Gene 2596 Human</a>
Synonyms:	Neuromodulin

#### Product images:



Dissociated cell culture from an e13 mouse brain showing GAP-43 (green staining) of a neuron. Note the presence of GAP-43 within the cell body, as well as within the growth cones at the end of each neurite.



[View online »](#)