


## Product datasheet for CH23008-100

### PLP1 Chicken Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	proteolipid protein 1
Database Link:	<a href="#">Entrez Gene 5354 Human</a>
Synonyms:	Myelin proteolipid protein, Lipophilin, PLP

#### Product images:



[View online »](#)

