


Product datasheet for CH23010-100

Vimentin (VIM) Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	vimentin
Database Link:	Entrez Gene 7431 Human
Synonyms:	VIM


Product images:



cultured dissociated cells from an e13 mouse brain), vimentin (green staining) can be seen forming the cytoskeleton of immature astrocyte progenitor cells. DAPI (blue staining) allows visualization of nuclei.



[View online »](#)



cultured dissociated cells from an e13 mouse brain), vimentin (green staining) can be seen forming the cytoskeleton of immature astrocyte progenitor cells. DAPI (blue staining) allows visualization of nuclei.