


## Product datasheet for CH23011-200

### GFAP Chicken Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse
Host:	Chicken
Isotype:	Ig
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens glial fibrillary acidic protein (GFAP), transcript variant 2
Database Link:	<a href="#">Entrez Gene 14580 Mouse</a> <a href="#">Entrez Gene 2670 Human</a>
Synonyms:	Glial Fibrillary Acidic Protein
Protein Families:	ES Cell Differentiation/IPS

#### Product images:



[View online »](#)

