


Product datasheet for CH23017-100

Coronin 1a (CORO1A) Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens coronin 1A (CORO1A), transcript variant 2
Database Link:	Entrez Gene 12721 Mouse Entrez Gene 11151 Human
Synonyms:	Coronin-like protein p57, Coronin-like protein, CORO1A, CORO1, Clipin-A, TACO

Product images:



[View online »](#)

