


Product datasheet for CH23019-200

Acpp Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Mus musculus acid phosphatase, prostate (Acpp), transcript variant 2
Database Link:	Entrez Gene 56318 Mouse
Synonyms:	Prostate acid phosphatase, PAP, ACP3, PSAP


Product images:



Tissue sections through a 10%-formalin fixed adult mouse. Sections were examined for PAP positive dorsal root ganglion sensory neurons. Secondary antibodies (fluorescein-labeled goat anti-chicken IgY, Neuromics Cat. #[CH23101]) were used at a concentration of 1:500.



[View online »](#)



Tissue sections through a 10%-formalin fixed adult mouse. Sections were examined for PAP positive dorsal root ganglion sensory neurons. Secondary antibodies (fluorescein-labeled goat anti-chicken IgY, Neuromics Cat. #[CH23101]) were used at a concentration of 1:500.