


## Product datasheet for CH23024-100

### Metabotropic Glutamate Receptor 2 (GRM2) Chicken Polyclonal Antibody

#### Product data:

|                |   |
|----------------|---|
| Product Type:  | Primary Antibodies                        |
| Applications:  | IF, IHC                                   |
| Reactivity:    | Human, Mouse                              |
| Host:          | Chicken                                   |
| Clonality:     | Polyclonal                                |
| Formulation:   | State: Aff - Purified                     |
| Gene Name:     | glutamate metabotropic receptor 2         |
| Database Link: | <a href="#">Entrez Gene 2912 Human</a>    |
| Synonyms:      | Metabotropic glutamate receptor 2, GPRC1B |

#### Product images:



Adult mouse cochlear ganglion stained with mGluR3 (green) and Tuj-1 (red).



[View online »](#)