


## Product datasheet for **CH23025-100**

### Metabotropic Glutamate Receptor 3 (GRM3) Chicken Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	glutamate metabotropic receptor 3
Database Link:	<a href="#">Entrez Gene 2913 Human</a>
Synonyms:	Metabotropic glutamate receptor 3, GPRC1C, MGLUR3

#### Product images:



Adult mouse cochlear ganglion stained with mGluR3 (green), Tuj-1 (red) and DAPI (blue).



[View online »](#)