


## Product datasheet for CH23026-100

### Grm5 Chicken Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Human, Mouse
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	glutamate receptor, metabotropic 5
Database Link:	<a href="#">Entrez Gene 108071 Mouse</a>
Synonyms:	Metabotropic glutamate receptor 5, GPRC1E, MGLUR5

#### Product images:



Adult mouse cochlear ganglion stained with mGluR5 (green) and Tuj-1 (red).



[View online »](#)