


Product datasheet for CH23027-200

Prion protein PrP (PRNP) Chicken Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF
Reactivity:	Human, Mouse, Rat
Host:	Chicken
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	prion protein
Database Link:	Entrez Gene 5621 Human
Synonyms:	Major prion protein, PrP27-30, PrP33-35C, ASCR, PRNP, PRIP


Product images:



Prion Protein staining of cultured fibroblasts cells transiently transfected with a plasmid containing the prion protein cDNA.



[View online »](#)



Prion Protein staining of cultured fibroblasts cells transiently transfected with a plasmid containing the prion protein cDNA.