

Product datasheet for **CL064**

MHC Class I Qa.m9 Mouse Monoclonal Antibody [Clone ID: 5034-54.8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	5034-54.8
Applications:	CT
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody reacts to Qa.m9.
Formulation:	State: Ascites State: Lyophilised Ascites
Conjugation:	Unconjugated
Database Link:	P14431



[View online »](#)