

## Product datasheet for CP507998

### Tissue Protein Lysate, Pancreas

#### Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Cancer
Diagnosis Category:	Carcinoid/Neuroendocrine Tumors
Tissue:	Pancreas
Frozen Sample ID:	FR000198E0
Case ID:	CI0000005744
Prep:	Not applic
Vial ID:	FR000198E0
Tissue of Origin:	Pancreas
Site of Finding:	Pancreas
Appearance:	Tumor
Sample Diagnosis (from Pathology Verification):	Tumor of pancreas, Islet cell
Normal:	0%
Lesional:	0%
Tumor:	80%
Tumor Hypo/Acellular Stroma:	20%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Malignant pancreatic endocrine cell tumor
CASE Diagnosis (from medical center path report):	Tumor of pancreas, neuroendocrine
Age:	60



[View online »](#)

<b>Gender:</b>	Female
<b>Race:</b>	Not reported
<b>Tumor Grade:</b>	Not Reported
<b>Minimum Stage Grouping:</b>	Not Reported
<b>T (Extent of Primary Tumor):</b>	Not Reported
<b>N (Lymph node metastasis):</b>	Not Reported
<b>M (Distant metastasis):</b>	Not Reported
<b>Diagnostic Test Results:</b>	Synaptophysin ~ SYN ! by IHC,Positive, diffuse. Chromogranin A ~ Parathyroid secretory protein 1 ! by IHC,Positive, focal
<b>Storage:</b>	Store the tissue lysate at -80°C.
<b>Stability:</b>	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.