

Product datasheet for CP508155

Tissue Protein Lysate, Colon

Product data:

| | |
|---|-------------------------|
| Product Type: | Tissue Protein Lysates |
| Disease State: | Cancer |
| Diagnosis Category: | Colorectal Cancer |
| Tissue: | Colon |
| Frozen Sample ID: | FR0001B2E3 |
| Case ID: | CI0000005822 |
| Prep: | Not applic |
| Vial ID: | FR0001B2E3 |
| Tissue of Origin: | Colon |
| Site of Finding: | Colon |
| Appearance: | Tumor |
| Sample Diagnosis (from Pathology Verification): | Adenocarcinoma of colon |
| Normal: | 0% |
| Lesional: | 0% |
| Tumor: | 70% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 25% |
| Necrosis: | 5% |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of colon |
| Age: | 71 |
| Gender: | Female |
| Race: | Not reported |


[View online »](#)

| | |
|-------------------------------------|--|
| Tumor Grade: | AJCC G2: Moderately differentiated |
| Minimum Stage Grouping: | I |
| T (Extent of Primary Tumor): | T2 |
| N (Lymph node metastasis): | N0 |
| M (Distant metastasis): | MX |
| Storage: | Store the tissue lysate at -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |