

## Product datasheet for CP511610

### Tissue Protein Lysate, Prostate

#### Product data:

|   |  |
|---|--|
| Product Type:                                     | Tissue Protein Lysates   |
| Disease State:                                    | Normal   |
| Diagnosis Category:                               | Normal   |
| Tissue:   | Prostate   |
| Frozen Sample ID:                                 | FR5B3390FC   |
| Case ID:  | CI0000006771   |
| Prep:   | Not applic   |
| Vial ID:  | FR5B3390FC   |
| Tissue of Origin:                                 | Prostate   |
| Site of Finding:                                  | Prostate   |
| Appearance:                                       | Normal   |
| Sample Diagnosis (from Pathology Verification):   | Within normal limits   |
| Normal:   | 100%   |
| Lesional:   | 0%   |
| Tumor:  | 0%   |
| Tumor Hypo/Acellular Stroma:                      | 0%   |
| Tumor Hypercellular Stroma:                       | 0%   |
| Necrosis:   | 0%   |
| Pathology Verification notes from H&E review:     | Non Tumor Structures: 50% Epithelium, 50% Fibromuscular stroma |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of prostate                                     |
| Age:  | 63   |



[View online »](#)

|                                     |  |
|-------------------------------------|--|
| <b>Gender:</b>                      | Male   |
| <b>Race:</b>                        | White or Caucasian   |
| <b>Tumor Grade:</b>                 | Not applicable   |
| <b>Minimum Stage Grouping:</b>      | Not applicable   |
| <b>T (Extent of Primary Tumor):</b> | Not applicable   |
| <b>N (Lymph node metastasis):</b>   | Not applicable   |
| <b>M (Distant metastasis):</b>      | Not applicable   |
| <b>Storage:</b>                     | Store the tissue lysate at -80°C.  |
| <b>Stability:</b>                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |