

Product datasheet for CP525908

Tissue Protein Lysate, Spleen

Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Spleen
Frozen Sample ID:	FR65EE21DA
Case ID:	CI0000009837
Prep:	Not applic
Vial ID:	FR65EE21DA
Tissue of Origin:	Spleen
Site of Finding:	Spleen
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Non Tumor Structures: 60% Red pulp, 40% White pulp
CASE Diagnosis (from medical center path report):	Lymphoma, large B-cell, diffuse
Age:	59



[View online »](#)

Gender:	Male
Race:	White or Caucasian
Tumor Grade:	Not applicable
Minimum Stage Grouping:	Not applicable
T (Extent of Primary Tumor):	Not applicable
N (Lymph node metastasis):	Not applicable
M (Distant metastasis):	Not applicable
Diagnostic Test Results:	CD3 ~ Leu4 ~ UCHT1 ! by IHC,Positive. Kappa light chain ! by FC,Positive, monoclonal. CD19 ~ B4 ~ Leu12 ! by IHC,Positive. CD22 ~ Leu14 ! by IHC,Positive. CD5 ~ Leu1 ! by IHC,Negative. CD10 ~ CALLA ~ Common ALL antigen ! by IHC,Negative. C
Storage:	Store the tissue lysate at -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.