

## Product datasheet for CP565393

### Tissue Protein Lysate, Ovary

#### Product data:

|   |                                 |
|---|---------------------------------|
| Product Type:                                     | Tissue Protein Lysates          |
| Disease State:                                    | Cancer                          |
| Diagnosis Category:                               | Ovarian Cancer                  |
| Tissue:   | Ovary                           |
| Frozen Sample ID:                                 | FR00002C16                      |
| Case ID:  | CI0000000069                    |
| Prep:   | 2P                              |
| Vial ID:  | PL00000779                      |
| Tissue of Origin:                                 | Ovary                           |
| Site of Finding:                                  | Ovary                           |
| Appearance:                                       | Tumor                           |
| Sample Diagnosis (from Pathology Verification):   | Adenocarcinoma of ovary, serous |
| Normal:   | 0%                              |
| Lesional:   | 0%                              |
| Tumor:  | 75%                             |
| Tumor Hypo/Acellular Stroma:                      | 0%                              |
| Tumor Hypercellular Stroma:                       | 15%                             |
| Necrosis:   | 10%                             |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of ovary, serous |
| Age:  | 51                              |
| Gender:   | Female                          |
| Race:   | White or Caucasian              |



[View online »](#)

|                                     |  |
|-------------------------------------|--|
| <b>Tumor Grade:</b>                 | FIGO G2: Moderately differentiated   |
| <b>Minimum Stage Grouping:</b>      | III  |
| <b>T (Extent of Primary Tumor):</b> | T3   |
| <b>N (Lymph node metastasis):</b>   | NX   |
| <b>M (Distant metastasis):</b>      | MX   |
| <b>Storage:</b>                     | Store the tissue lysate at -80°C.  |
| <b>Stability:</b>                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |