

## Product datasheet for CP565470

### Tissue Protein Lysate, Kidney

#### Product data:

|   |  |
|---|--|
| Product Type:                                     | Tissue Protein Lysates                       |
| Disease State:                                    | Normal                                       |
| Diagnosis Category:                               | Normal                                       |
| Tissue:   | Kidney                                       |
| Frozen Sample ID:                                 | FR0001922E                                   |
| Case ID:  | CI0000005951                                 |
| Prep:   | 4P   |
| Vial ID:  | PL00000549                                   |
| Tissue of Origin:                                 | Kidney                                       |
| Site of Finding:                                  | Kidney                                       |
| Appearance:                                       | Normal                                       |
| Sample Diagnosis (from Pathology Verification):   | Within normal limits                         |
| Normal:   | 100%   |
| Lesional:   | 0%   |
| Tumor:  | 0%   |
| Tumor Hypo/Acellular Stroma:                      | 0%   |
| Tumor Hypercellular Stroma:                       | 0%   |
| Necrosis:   | 0%   |
| Pathology Verification notes from H&E review:     | 100% cortex; focal mild chronic inflammation |
| CASE Diagnosis (from medical center path report): | Carcinoma of kidney, renal cell, clear cell  |
| Age:  | 74   |



[View online »](#)

|                                     |  |
|-------------------------------------|--|
| <b>Gender:</b>                      | Male   |
| <b>Race:</b>                        | White or Caucasian   |
| <b>Tumor Grade:</b>                 | Not applicable   |
| <b>Minimum Stage Grouping:</b>      | Not applicable   |
| <b>T (Extent of Primary Tumor):</b> | Not applicable   |
| <b>N (Lymph node metastasis):</b>   | Not applicable   |
| <b>M (Distant metastasis):</b>      | Not applicable   |
| <b>Storage:</b>                     | Store the tissue lysate at -80°C.  |
| <b>Stability:</b>                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |