

Product datasheet for CP565526

Tissue Protein Lysate, Esophagus

Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Esophagus
Frozen Sample ID:	FR00026911
Case ID:	CI0000008013
Prep:	2P
Vial ID:	PL000028EC
Tissue of Origin:	Esophagus
Site of Finding:	Esophagus
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	5% mucosa, 20% submucosa, 75% muscle
CASE Diagnosis (from medical center path report):	Adenocarcinoma of stomach, signet ring cell
Age:	71



[View online »](#)

Gender:	Male
Race:	White or Caucasian
Tumor Grade:	Not applicable
Minimum Stage Grouping:	Not applicable
T (Extent of Primary Tumor):	Not applicable
N (Lymph node metastasis):	Not applicable
M (Distant metastasis):	Not applicable
Storage:	Store the tissue lysate at -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.