

Product datasheet for CP565539

Tissue Protein Lysate, Lung

Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Other Disease
Diagnosis Category:	Pneumonia
Tissue:	Lung
Frozen Sample ID:	FR000270F8
Case ID:	CI0000008283
Prep:	2P
Vial ID:	PL00002262
Tissue of Origin:	Lung
Site of Finding:	Lung
Appearance:	Lesional
Sample Diagnosis (from Pathology Verification):	Bronchiolitis obliterans organizing pneumonia
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	no residual normal architecture; 2+ MAC, 2+ LYP
CASE Diagnosis (from medical center path report):	Adenocarcinoma of lung
Age:	75



[View online »](#)

Gender:	Female
Race:	White or Caucasian
Tumor Grade:	Not applicable
Minimum Stage Grouping:	Not applicable
T (Extent of Primary Tumor):	Not applicable
N (Lymph node metastasis):	Not applicable
M (Distant metastasis):	Not applicable
Storage:	Store the tissue lysate at -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.