

Product datasheet for CP565609

Tissue Protein Lysate, Stomach

Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Cancer
Diagnosis Category:	Stomach Cancer
Tissue:	Stomach
Frozen Sample ID:	FR0002D6A0
Case ID:	CI0000010335
Prep:	1P
Vial ID:	PL00000266
Tissue of Origin:	Stomach
Site of Finding:	Stomach
Appearance:	Tumor
Sample Diagnosis (from Pathology Verification):	Adenocarcinoma of stomach, signet ring cell
Normal:	0%
Lesional:	0%
Tumor:	75%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	25%
Necrosis:	0%
Pathology Verification notes from H&E review:	Tumor Stroma (Cellular): Desmoplastic reaction, Inflammatory cells
CASE Diagnosis (from medical center path report):	Adenocarcinoma of stomach, signet ring cell
Age:	60



[View online »](#)

Gender:	Male
Race:	White or Caucasian
Tumor Grade:	AJCC G3: Poorly differentiated
Minimum Stage Grouping:	IV
T (Extent of Primary Tumor):	T2
N (Lymph node metastasis):	N2
M (Distant metastasis):	M1
Diagnostic Test Results:	Alcian blue stain,Positive. Mucicarmine,Positive. Cytokeratin 7 ~ CK7 ! by IHC,Positive. Cytokeratin 20 ~ CK20 ! by IHC,Positive, focal. Cytokeratin CAM5.2 ~ CAM5.2 ! by IHC,Positive. Cytokeratin AE1/AE3 ~ AE1/AE3 ! by IHC,Positive
Storage:	Store the tissue lysate at -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.