

Product datasheet for CP565610

Tissue Protein Lysate, Stomach

Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Stomach
Frozen Sample ID:	FR0002D6AC
Case ID:	CI0000010335
Prep:	2P
Vial ID:	PL00000283
Tissue of Origin:	Stomach
Site of Finding:	Stomach
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Non Tumor Structures: 25% Mucosa, 40% Submucosa, 35% Muscle
CASE Diagnosis (from medical center path report):	Adenocarcinoma of stomach, signet ring cell
Age:	60



[View online »](#)

Gender:	Male
Race:	White or Caucasian
Tumor Grade:	Not applicable
Minimum Stage Grouping:	Not applicable
T (Extent of Primary Tumor):	Not applicable
N (Lymph node metastasis):	Not applicable
M (Distant metastasis):	Not applicable
Diagnostic Test Results:	Alcian blue stain,Positive. Mucicarmine,Positive. Cytokeratin 7 ~ CK7 ! by IHC,Positive. Cytokeratin 20 ~ CK20 ! by IHC,Positive, focal. Cytokeratin CAM5.2 ~ CAM5.2 ! by IHC,Positive. Cytokeratin AE1/AE3 ~ AE1/AE3 ! by IHC,Positive
Storage:	Store the tissue lysate at -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.