

Product datasheet for CP565649

Tissue Protein Lysate, Colon: right

Product data:

| | |
|---|---|
| Product Type: | Tissue Protein Lysates |
| Disease State: | Normal |
| Diagnosis Category: | Normal |
| Tissue: | Colon |
| Frozen Sample ID: | FR00037A07 |
| Case ID: | CI0000013631 |
| Prep: | 1P |
| Vial ID: | PL00002321 |
| Tissue of Origin: | Colon: right |
| Site of Finding: | Colon: right |
| Appearance: | Normal |
| Sample Diagnosis (from Pathology Verification): | Within normal limits |
| Normal: | 100% |
| Lesional: | 0% |
| Tumor: | 0% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 0% |
| Necrosis: | 0% |
| Pathology Verification notes from H&E review: | Non Tumor Structures: 30% Mucosa, 40% Submucosa, 30% Muscle |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of colon |
| Age: | 74 |



[View online »](#)

| | |
|-------------------------------------|--|
| Gender: | Female |
| Race: | Black or African American |
| Tumor Grade: | Not applicable |
| Minimum Stage Grouping: | Not applicable |
| T (Extent of Primary Tumor): | Not applicable |
| N (Lymph node metastasis): | Not applicable |
| M (Distant metastasis): | Not applicable |
| Storage: | Store the tissue lysate at -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |